



rubitec ag


VARUMILL Solid calibrating for the chemical and pharmaceutical industries. No off-the-peg item but custom-tailored to your individual need

VARUMILL combines ideally the proven conical and cylindrical arrangements of crushing screens and rotor with the technical advantages of the rubitec system.

- CIP-compatible even with inbuilt screen
- design conforming to GMP
- energy supply from technical room according to GMP
- enclosed and cooled drive for clean room operation
- free choice of inlet and outlet branches
- no screws in production cavity
- use in Ex -zone



Optionally: 1 or 2 additional rotors for poorly flowing powder



- additional rotor to remove powder between sieve and housing
- additional rotor to prevent powder bridge in in feed hopper of pneumatic conveyer

CH-4431 Bennwil, Hauptstrasse 8
Tel +41 -(0)61-956 91 91
Fax +41 -(0)61-956 91 90
info@rubitec.ch
www.rubitec.ch

rubitec ag
chemical-pharmaceutical-equipment



- ▮ Screens and rotors for different products: dry or moist
- ▮ Screens changed without tools
- ▮ Gas-flushed shaft seal



rubitec ag

Our engineering department will be pleased to assist you in integrating **VARUMILL** units into your overall plant planning.

And:

For emergencies or routine maintenance, our service organization is on call with original rubitec spare parts.



rubitec ag



Plant and mechanical engineering

Explosion Isolation

Solid Calibration

Lifting/Tilting Equipment

- ▮ Use of conical and cylindrical screens with appropriate rotors



CH-4431 Bennwil, Hauptstrasse 8
Tel +41 -(0)61-956 91 91
Fax +41 -(0)61-956 91 90
info@rubitec.ch
www.rubitec.ch

rubitec ag
chemical-pharmaceutical-equipment